

r26m PT

Diagnostic Tympanometer



Compact. Ergonomic. User friendly.

The Resonance® R26M Middle Ear Analyzer provides a flexible clinical device for tympanometry acoustic reflex measurement. Compact dimensions, excellent ergonomics and a luminous 7" LCD display, the R26M has been engineered to respond to the needs of ENT and Audiology Specialists. Design is user-friendly and allows both device and accessories to be safely transported in a durable luggage-grade padded carrying case, available on request.

Safe. Flexible. Customizable.

The R26M offers several options for managing and archiving PDF test-report files. Test-report files, generated in PDF, can be initially stored in a USB-stick, then archived into external units or printed thanks to the built-in thermal printer. Reports can be easily edited and customized with additional text and logos. The Management Data Suite (MDS) software allows patient data and test results to be transferred to the PC and to Noah®



Fast. Versatile. Qualitative.

The R26M is built upon 30 years of experience and know-how gained in the audiology industry by Resonance. The R26M configuration offers a complete range of diagnostic tympanometry tests. The R26M PT configuration includes Multifrequency Tympanometry. Both configurations aim to simplify daily operations and maximize efficiency.

Comfortable. Elegant. Reliable.

The design of accessories maximizes patient comfort for the best and easiest handling of the device. Optional accessories increase flexibility and performance. Resonance R&D collaborates with universities and ENT specialists to identify and prioritize features that are most relevant and important to them. High-quality components ensure device reliability as well as long-term operation.



Innovative thinking
has arrived.



Exclusive U.S. Distributor:
EAR TECHNOLOGY CORPORATION
P.O. Box 1516
Johnson City, TN 37605 USA
1-800-327-8547
resonance-eartech.com | resonance@eartech.com

GENERAL SPECIFICATIONS

TEST TYPES

- Tympanometry: Automatic Tympanometry, Acoustic Reflex, Reflex Decay, Multifrequency Tympanometry

DIMENSIONS AND WEIGHT

- L x W x H: 14.6" x 11.4" x 7.1"
- Net Weight: 5.5 lbs.
- Body material: Bayblend® self-extinguishing

DISPLAY

- 7" LCD color display

USER INTERFACE

- Multilingual

PRINTER

- Built-in fast thermal printer with paper width: 4.4" supplied as standard part

REPORTS

- Printed on thermal paper
- PDF report created directly from the device and stored on USB pen drive with possibility to add patient data and test comments via the USB Keyboard (optional)
- Data transfer to PC using MDS* software

"CHILDREN" FEATURE

- To help keep the child's attention diverted from the testing procedure, a series of animated images appears on the color display

COMMUNICATION PORT

- Nr.1 USB host type A
- Nr.1 USB slave type B

PC INTERFACE

- MDS* software with reporting and printing features
- Compatible with Windows® 7, 8, 8.1, 10
- Data transfer via USB cable

POWER

POWER SUPPLY

- 110-240V AC 50/60Hz 40VA
- Fuses: 2 x T 1 A L 250V

CONSUMPTION

- Max current: 0.15A
- Power consumption: 40VA

QUALITY SYSTEM

Manufactured, designed and developed under an ISO 13485, ISO 9001 certified quality system. Medical CE mark and FDA approval.

COMPLIANCE/REGULATORY STANDARDS

- Certifying body: IMQ S.p.A. Istituto Italiano Del Marchio Di Qualita
- MDD 93/42/EEC and update 2007/47/CEE: Class IIa (as referred to in Annex IX, rule 10)
- Safety: IEC 60601-1, 3rd and 2nd edition, Class 1 Type B
- EMC: IEC 60601-1-2
- Impedance: IEC 60645-5/ANSI S.3.39 Type 2

OPTIONAL

- External USB keyboard
- Durable luggage-grade padded carrying case
- Foot switch
- IP30 Insert earphones for contra use

TYMPANOMETRY SPECS

PROBE TONE

- 226 Hz for Admittance (Y) curve tympanometry
- 678, 800 and 1000 Hz for multifrequency tympanometry with added "B" and "G" curve

INTENSITY

- 226 Hz: 85dB SPL \pm 2dB
- 678, 800 and 1000 Hz: 75dB SPL \pm 2dB
- Frequency Accuracy: \pm 0.5%
- Harmonic distortion: < 1%

ADMITTANCE MEASUREMENTS

- Compliance range at 226Hz: 0.05 up to 7ml
- Sensitivity scale: Autoscale to appropriate range; available scale at 226Hz: 1.5, 2, 5ml or 7ml
- 678, 800 and 1000 Hz: 0 to +25 mmho
- Sensitivity scale HF: Autoscale to appropriate range, available scales at 678, 800 and 1000Hz: 5, 10, 15, 20, 25 mmho

AIR PRESSURE

- Control: Automatic
- Range: from +200 up to -400 daPa, both adjustable in 50 daPa steps
- Sweep rate: 50, 100, 200, 300 daPa/sec and automatic
- Indicator: Measured value is displayed
- Safety limitations: -800 up to +600 daPa
- Pressure accuracy: \pm 10 daPa or \pm 10%

ACOUSTIC REFLEX TESTS

- Reflex test method: Threshold, Automatic, Manual
- Stimulus duration: 0.5, 1 or 2 sec.
- User selectable protocols for all test methods
- Ipsi or contralateral stimulation for all reflex test
- Automated peak search functions available for all test methods
- Manual Reflex: manual control of all stimuli
- Reflex Decay: Threshold, Automatic, Manual
- Stimulus duration: 10 or 20 sec.
- Stimulus duration: 1 sec. fixed

FREQUENCIES AND INTENSITY RANGES

- Ipsilateral:
 - Pure tone level range (dB HL) from 50 to 110 (all versions)
 - Noise level range (dB SPL) from 50 to 100 (all versions)
- Frequency: 500, 1000, 2000, 3000, 4000 Hz, BBN, HP or LP noise
- Accuracy: \pm 1%
- Harmonic distortion (THD): < 3%
- Contralateral:
 - Pure tone level range (dB HL) from 50 to 120 (all version)
 - Noise level range (dB SPL) from 50 to 115 (all version)
- Frequency: 250, 500, 1000, 2000, 3000, 4000, 6000, 8000 Hz, BBN, HP or LP noise
- Accuracy: \pm 1%
- Calibration accuracy \pm 3 dB
- Level steps: 1, 2, 5 or 10 dB
- On/Off ratio: 70 dB minimum

STANDARD ACCESSORIES

- Headband Probe 226 Hz or Headband Probe HF with DD45 for contra
- Insert for contra
- Kit of assorted ear tips; 2 shapes available
- Calibration cavity with probe housing
- Built-in fast thermal printer
- Thermal paper roll
- Device dust cover
- USB pen drive
- MDS* software
- Power supply cable (110 – 220V)

* Resonance Management Data Suite

Exclusive U.S. Distributor:

EAR TECHNOLOGY CORPORATION

P.O. Box 1516
Johnson City, TN 37605 USA
1-800-327-8547

resonance-eartech.com | resonance@eartech.com

